HYALEXO CROSS

Hyaluronate de Sodium réticulé

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DESCRIPTION

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INSTRUCTIONS

DESCRIPTION

HYALEXO CROSS is a sterile, apyrogenic, resorbable viscoelastic solution made with cross-linked hyaluronic acid created from non-animal sourced bacterial bio-fermentation and administered via intra-articular injection. HYALEXO CROSS is a colourless gel contained in a pre-filled, graduated, disposable syringe with a Luer-Lok adapter, which contains 3 ml of solution.

COMPOSITION

1 ml of solution contains: 20mg of cross-linked sodium hyaluronate, sterile saline solution.

RECOMMENDED USE

Pain and reduced mobility due to traumatic or degenerative deformities of the knee joint and other synovial joints. Non-toxic to animals.

MECHANISM OF ACTION

Hyaluronic acid is naturally synthesised by synovial cells. In addition to being responsible for the viscoelastic properties of synovial fluid, it helps lubrication mechanisms under loading conditions, and has a protective function. In degenerative joint diseases, such as chronic primary osteoarthritis, synovial fluid has a reduced viscosity, decreasing its ability to lubricate the joint and absorb shocks. Viscosupplementation with hyaluronic acid improves mobility and joint function, thereby reducing mechanical joint overload and, consequently, painful symptoms.

CONTRAINDICATIONS

HYALEXO CROSS is not recommended for patients with acknowledged hypersensitivity to the product ingredients. HYALEXO CROSS must not be injected into the joint in the event of venous or lymphatic stasis.

MEDICAL WARNING

HYALEXO CROSS is packaged to be used only once, on one patient during one session. HYALEXO CROSS is a sterile, disposable product and must not be resterilised. Use HYALEXO CROSS straight after opening the blister and dispose of the syringe immediately after use, even if the entire content has not been used. Do not use the product after the expiry date printed on the package. As there are no clinical data available on the use of hyaluronic acid in children, during pregnancy, or breastfeeding, treatments with HYALEXO CROSS are not recommended in these instances. Not to be used on people below the age of eighteen. Do not inject endovascularly, aspirating with the syringe prior to the injection. Do not inject in extra-articular sites or in the thickness of adjacent tissue. HYALEXO CROSS must not be used in the event of significant intra-articular effusion prior to the injection.

Keep the product out of reach from children. HYALEXO CROSS must not be used for purposes other than those described in this leaflet.

SAFETY PRECAUTIONS

HYALEXO CROSS is a colourless transparent gel. If the content of the syringe appears cloudy or contains particles, do not use the product and promptly notify Apharm s.r.l. at the address below.

Do not use the device if the sterile syringe or the blister containing it appear to be damaged. Dispose of the syringes and needles after using the medical device in compliance with the approved medical practice.

Have the patient avoid demanding physical activities for the joint for at least a couple of days after the treatment. Do not inject the product in the event of an infected or severely swollen joints – for example in the event of rheumatoid arthritis or Bechterew's disease – or if the patient has an infection in the area of the treatment. Do not inject the product in the event of skin lesions, skin redness, or any type of dermatitis.

UNDESIRABLE EFFECTS

An HYALEXO CROSS joint infiltration may produce local secondary effects, such as pain, feeling of heat, redness, and swelling, mainly caused by the trauma of the needle. These symptoms can be alleviated by applying ice on the treated part for 5-10 minutes and by resting the joint.

IMPORTANT: the physician must ensure that the patient reports any type of adverse reaction, whether serious or not, even if it is not listed in this information leaflet and must promptly notify Apharm s.r.l. at the address below.

INCOMPATIBILITY: there is a known incompatibility between hyaluronic acid and quaternary ammonium salts, such as benzalkonium chloride. Therefore, it is important to avoid placing HYALEXO CROSS into contact with these substances.

DIRECTIONS

HYALEXO CROSS must be injected into the joint cavity using a spinal needle ranging between 18 and 21G. Use a needle of a suitable length for treating the joint. The intra-articular infiltration must be performed solely by qualified medical personnel and authorised in compliance with the laws in force. Take all the normal precautions for intra-articular injections. This product must be used with caution in the event of lymphatic or venous stasis in the leg to be treated. In the event of joint effusion, we recommend reducing the effusion through aspiration, rest, and application of an ice pack. In this instance, we recommend starting the HYALEXO CROSS treatment 2 or 3 days later. Inject the solution contained in HYALEXO CROSS making sure that it is at room temperature. To remove the syringe from its protective blister, keep it close to the body without touching the stem of the plunger. Unscrew the cap of the syringe and place the needle onto the Luer-Lok adapter by screwing on the base clockwise until it locks into place. For a safe use of HYALEXO CROSS, it is important to assemble the needle and syringe correctly, as an incorrect assembly may cause the needle to separate from the syringe during the injection (see figures 1 and 2).

Remove the nozzle of the syringe with caution using aseptic techniques.

Remove the protective cap of the needle only just before making the injection.

Before injecting, expel the air from the syringe until a drop of the product appears on the tip of the needle.

Inject the product in perfectly intact skin, after removing any cosmetic residue and disinfecting it.

Inject the product only in the synovial space. If necessary, carry out this operation with a guiding instrument, e.g. fluoroscopic guidance, especially in the event of hip and shoulder treatments. Ionic or non-ionic contrast media can be used under fluoroscopic guidance.

Local anaesthetic infiltration is recommended prior to the infiltration, especially in the event of hip infiltration.

Fig. 1 Unscrew the cap of the Luer-Lok adapter


Fig. 2 Insert the needle

DOSAGE

Perform a single injection of HYALEXO CROSS in the affected joint. Do not exceed the quantity of one syringe (3 ml) per injection site and per session. The beneficial effect of a treatment cycle can last up to 6 months, according to the severity of the joint condition.

In the event of additional treatment cycles, we recommend a 6-month interval between treatments. The HYALEXO CROSS treatment affects only the joint with the injection site and does not produce a general systemic effect.

PACKAGE

One packaging of HYALEXO CROSS contains:

- One 3 ml pre-filled, graduated glass syringe with a Luer-Lok adapter sealed in a sterile packaging.
- One set of adhesive labels with the lot number, one of which must be applied to the patient's clinical file, and the other to be delivered to the patient. These procedures are compulsory to ensure the traceability of the product.
- Directions complying with European regulations to be read carefully before using the device.

STORAGE INSTRUCTIONS

Read the expiry date of the product on the syringe, on the label on the blister, or on the packaging. Store in a dry place, away from light and at a temperature ranging between 5 °C and 25 °C. Do not freeze.

Note: the graduation on the syringe is not for measurement purposes but is only provided as a visual control over the injected material.

SYMBOLS ON THE PACKAGE

	Do not reuse, disposable product
	Carefully follow the instructions
	Expiry date
	Lot number
	Sterile for aseptic manufacturing
	Keep in a dry place
	Shield from sunlight
	Store at a temperature between 5 °C and 25 °C
	Manufacturer
	Date of Manufacture
	Product code
	Do not use if packaging is damaged
	CE trademark compliant with 93/42/EEC concerning medical devices. 0373 is the number of the Notified Body

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