140*297mm

***** OK DEF *****

HYALEXO EXPANSCIENCE

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Produit

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Fabricant

Product

Hyalexo Expanscience is a CE-marked medical device fulfilling the requirements of MDD 93/42/EEC.

Composition

Principal ingredient: sodium hyaluronate (EP) 1%;
Other ingredients: sodium chloride, dibasic sodium phosphate hydrate, monobasic sodium phosphate, water for injection.

Description

Hyalexo Expanscience is a sterile, clear colorless viscoelastic solution supplied in a glass syringe. The solution having a volume of 2 ml contains 20 mg of sodium hyaluronate of very high molecular weight fraction of sodium hyaluronate (2500-3000 kDa) obtained from bacterial fermentation. Hyaluronic acid is a polysaccharide naturally present in several human tissues such as cartilage.

It is the major constituent of synovial fluid and provides its viscosity and elasticity. In traumatic and degenerative joint diseases, the amount of hyaluronic acid is reduced. The consequent loss of viscosity by the synovial fluid may be the cause of painful symptoms. Intra-articular administration of sodium hyaluronate may restore the viscoelastic properties of synovial fluid, attenuate pain and improve the mobility of the joint. Product is almost biodegraded within 4 weeks after implantation.

Indications

Hyalexo Expanscience is indicated for the symptomatic relief of pain in patients suffering from osteoarthritis (OA) or degenerative joint diseases causing an alteration of the functionality of the synovial fluid.

Dosage and administration

Each syringe is a single dose preparation. Inject one syringe of Hyalexo Expanscience using a sterile needle (18 to 22G) into the affected joint once a week for three weeks using a standard technique. Hyalexo Expanscience should be administered only by qualified physicians and following strict aseptic conditions. If joint effusion is present, it should be aspirated by arthrocentesis prior to the injection.

Contraindication

- 1) Do not use to treat patients with known hypersensitivity to the principal ingredient or other components of the product.
- 2) Do not use in the presence of skin infections or skin diseases in the area of the injection site.
- 3) Do not use if the knee joint is infected or severely inflamed.

Adverse reactions

- 1) Local reactions at the injection site may occur such as discomfort, transient pain (i.e. arthralgia, myalgia,...) and transient increase of joint inflammation right after or within 48 hours after injection, joint swelling, joint effusion, redness, warm feeling, bleeding and ecchymosis. These side effects are generally of mild to moderate intensity and resolve spontaneously within a few days. Use aseptic technique when handling the syringe and administering the injection to avoid any risk of septic arthritis (joint infections).
- 2) Hypersensitivity to one of the components: Local allergic reactions such as rash (redness, skin eruption), local urticaria (hives) and itching may occur. These side effects are rare, of mild to moderate intensity and most often resolve spontaneously within a few days.

HYALEXO EXPANSCIENCE

(Pre-filled Syringe)

In very rare cases, systemic allergic reaction may occur with angioedema (i.e. facial edema, eyelid edema, bumps of tongue...), facial flushing, dyspnea (shortness breath), generalized oedema or urticaria. If such adverse reaction occurs, Hyalexo Expanscience should be discontinued and appropriate measures should be taken. In some cases, injection itself may cause stress-induced reaction such as nausea, vomiting and syncope.

Precautions

Hyalexo Expanscience should only be administered via intraarticular injection. Safety and effectiveness of Hyalexo Expanscience in joints other than the knee and shoulder have not been investigated. Co-administration of Hyalexo Expanscience with any other intra-articular treatment has not been established and should be avoided to prevent possible interaction. Do not use concomitantly with disinfectants containing quaternary ammonium salts such as benzalkonium chloride because hyaluronic acid may precipitate in their presence. Each syringe must be disposed of after its use. Do not use after expiration date or if the blister package is opened or damaged.

Pregnancy and Nursing mothers

- 1) Pregnancy: Animal reproduction studies have revealed no evidence of impaired fertility. However, Safety and efficacy have not been determined during pregnancy. Women who are pregnant or considering pregnancy should be administered only if the therapeutic benefit justifies the potential risks.
- 2) Nursing Mothers: It is not known if Hyalexo Expanscience is excreted in human milk. The safety and effectiveness of Hyalexo Expanscience have not been established in lactating women.

Pediatric use

The safety and effectiveness of Hyalexo Expanscience have not been tested in children.

Storage

Store at a temperature of 1°C – 30°C. Protect from direct sunlight and freezing.

Shelf life

Expiration date is indicated on the package and refers to the product stored in its original packaging according to the above storage conditions.

Packaging Unit

One cardboard package containing one 3mL pre-filled syringe contains 2ml of solution closed with a rubber stopper. The content of the syringe is sterile, and has been sterilized by aseptic processing techniques. The pre-filled syringe is packaged in a sealed blister.

Manufacturer

YooYoung Pharmaceutical
33, Yongso 2-gil, Gwanghyewon-myeon, Jincheon-gun, Chungcheongbuk-do, Korea
Tel: +82-2-6207-6114 / Fax: +82-2-6202-7067

Authorized Representative in the European Community

Obelis s.a.
Bd. Général Wahis, 53, B-1030 Brussels, Belgium
Tel: +32.2.732.59.54
mail@obelis.net

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Consult instructions for use

Single use only, do not re - Use

Temperature limitation (Store at a temperature of 1°C – 30°C)

Caution

Keep away from sun light

Date of manufacture

Authorized representative in the EC

Manufacturer

"Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established"

Made in Korea

Rev.003

Revised date 2021. 10. 05

FRPR2362110

Y
PHARMACEUTICAL
YOOYOUNG

We care for your packaging

Manufactured by : Yoo Young Pharm, Co., Ltd,

33,Yongso 2-gil, Gwanghyewon-myeon,
Jincheon-gun, Chungcheongbuk-do, Korea

APPROVAL PROOF N°

DATE

SIGNATURE

COLORS USED

Black

ARTWORK INFORMATION

Customer

ID record product

Code No.

Product Name

Dimension etuis

140*297